

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

ACTN1 Antibody (Clone AT1D10)

单克隆抗体

产品货号: MA01010

第三版

免疫原:Recombinant human ACTN1 (1-249aa) purified from E. coli

亚型:Mouse IgG2a, k

克隆:Anti-human ACTN1 mAb, clone AT1D10, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human ACTN1 protein.

描述:ACTN1 also known as Actinin, alpha 1. ACTN1 functions as an actin-binding-protein. ACTN1 determines the motility of keratinocytes by regulating the organization of the actin cytoskeleton, focal adhesion, and hemidesmosome proteins complexes, thereby modulating cell speed, lamellipodial dynamics, and directed migration.

配方:Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度:1 mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C 或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF