

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

AGXT Anti-Human Antibody (Clone AT2T4)

单克隆抗体

产品货号: MA01022

第三版

**免疫原:**Recombinant Maspin (1-375aa) purified from E. coli

**亚型:**Mouse IgG2b heavy chain and kappa light chain

**克隆:**Anti-human Maspin mAb, clone AT2N6, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human Maspin protein.

**描述:**The antibody has been tested by ELISA, Western blot analysis and Immunofluorescence analysis to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results. Recommended dilution range for Western blot analysis is 1:1000 ~ 2000.

**配方:**Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

**浓度:**1 mg/ml

**保存:** +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

**应用:**ELISA, WB, ICC/IF