

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

AIMP1 Antibody (Clone AT6D10)

单克隆抗体

产品货号: MA01025

第三版

**免疫原:** Recombinant human AIMP1 (1-336aa) purified from E. coli

**亚型:** Mouse IgG2a, k

**克隆:** Anti-human AIMP1 mAb, clone AT6D10, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human AIMP1 protein.

**描述:** AIMP1, also known as p43, hasn't until recently been recognized for its prominent immunological functions. Together with other nonenzymatic mARS structural components AIMp2/38 and AIMp3/p18, it participates in the machinery responsible for cell-cycle control and tumor suppression.

**配方:** Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

**浓度:** 1 mg/ml

**保存:** +4 ° C 保存 (1-2 周). 长期保存在-20 ° C 或者-70 ° C. 避免反复冻融.

**应用:** ELISA, WB, ICC/IF