

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

AK3 Anti-Human Antibody (Clone SJB3-36)

单克隆抗体

产品货号: MA01028

第三版

免疫原: Recombinant human AK3 (adenylate kinase isozyme 3) purified from E. coli

亚型: Mouse IgG1 heavy chain and kappa light chain

克隆: Anti-human AK3 mAb, clone SJB3-36, is derived from hybridization of mouse SP2/O myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human AK3 protein.

描述: Adenylate kinase (AK; adenosine triphosphate-adenosine monophosphate [ATP-AMP] phosphotransferase, EC 2.7.4.3) is a ubiquitous monomeric enzyme involved energy metabolism of prokaryotic and eukaryotic cells.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

应用: ELISA, WB, ICC/IF