



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

ARF1 Anti-Human Antibody (Clone AT1B3)

单克隆抗体

产品货号: MA01094

第三版

免疫原: Recombinant human ARF1 (1-181aa) purified from E. coli

亚型: Mouse IgG1 heavy chain and light chain

克隆: AT1B3

描述: ADP-ribosylation factor 1 (ARF1) is a member of the human ARF gene family. The family member is essential and ubiquitous in eukaryotes, being involved in vesicular transport and functioning as an activator of phospholipase D and cholera toxin.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1 mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

应用: ELISA, WB, ICC/IF