

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

ATGL Anti-Human Antibody (Clone AT2G2)

单克隆抗体

产品货号: MA01101

第三版

免疫原: Recombinant human ATGL (30-504aa) purified from E. coli

亚型: Mouse IgG2b heavy chain and kappa light chain

克隆: Anti-human ATGL mAb, clone AT2G2, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice

描述: Adipose triglyceride lipase (ATGL) is a 504 amino acid protein that is highly expressed in mouse and human adipose tissue. ATGL catalyzes the initial step in triglyceride hydrolysis in adipocyte lipid droplets and has acylglycerol transacylase activity. Inhibition of ATGL markedly decreases total adipose acyl-hydrolase activity.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1 mg/ml

保存: +4 °C 保存 (1-2 周). 长期保存在 -20 °C 或者 -70 °C. 避免反复冻融.

应用: ELISA, WB, ICC/IF