



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

AURKB Anti-Human Antibody (Clone AT2B1)

单克隆抗体

产品货号: MA01102

第三版

免疫原: Recombinant human Aurora kinase B (1-344aa) purified from E. coli

亚型: Mouse IgG2b heavy chain and kappa light chain

克隆: AT2B1

描述: The aurora family (A, B and C) are serine threonine kinases and key regulators of chromosome segregation during mitosis. Aurora kinase B is a chromosomal passenger protein that regulates chromosome segregation and cytokinesis.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1 mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

应用: ELISA, WB, ICC/IF