



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

BAK1 Anti-Human Antibody (AT8B4)

单克隆抗体

产品货号: MA01106

第三版

**免疫原:** Recombinant human BAK1 (29-187aa) purified from E. coli

**亚型:** Mouse IgG2a heavy chain and kappa light chain

**克隆:** Anti-human BAK1 mAb, clone AT8B4, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice

**描述:** BAK1 encodes an receptor-like kinase (RLK) with a putative extracellular domain, a single transmembrane domain, an intracellular-juxtamembrane domain, and a kinase domain. BAK1 expression is associated with the progression of Prostate cancer (CaP).

**配方:** Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

**浓度:** 1 mg/ml

**保存:** +4 °C 保存 (1-2 周). 长期保存在 -20 °C 或者 -70 °C. 避免反复冻融.

**应用:** ELISA, WB, ICC/IF