
本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Bcl-1/Cyclin D1 (5D4) Mouse IgG MoAb Biotin

单克隆抗体

产品货号: MA01109

第三版

免疫原:recombinant Human PRAD1/cyclin D1 (E. coli)

亚型:IgG2a

克隆:5D4

描述:Bcl-1/Cyclin D1 belongs to the G1 cyclins and plays a key role in cell cycle regulation during the G1/S transition by cooperating with cyclin-dependent kinases (CDKs).

配方:Lyophilized product from 1% BSA in PBS containing 0.05%NaN3

浓度:100 ug/mL

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF