

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Bcl-2 Anti-Human Antibody (AT1B5)

单克隆抗体

产品货号: MA01110

第三版

免疫原: Recombinant human Bcl-2 (1-211aa) purified from E. coli

亚型: Mouse IgG2b heavy chain and kappa light chain

克隆: Anti-human Bcl2 mAb, clone AT1B5, is derived from hybridization of mouse FO myeloma cells with spleen cells from BALB/c mice

描述: Bcl-2 (B-cell lymphoma protein 2 alpha) is an anti-apoptotic protein located primarily in the outer mitochondrial membrane. It blocks the apoptotic death of some cells such as lymphocytes. BCL-2 is thought to regulate cell death by controlling the mitochondrial membrane permeability during apoptosis.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1 mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

应用: ELISA, WB, ICC/IF