

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Blys Antibody (Clone H4-C7)

单克隆抗体

产品货号: MA01117

第三版

**免疫原:** Recombinant human BLYS (134-285aa) purified from E. coli

**亚型:** Mouse IgG2a heavy chain and kappa light chain.

**克隆:** Anti-human BLYS mAb, clone H4-C7, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice.

**描述:** B lymphocyte stimulator (BLYS), a member of the tumor necrosis factor (TNF) superfamily of cytokines, induces B-cell proliferation and immunoglobulin secretion and is a key regulator of peripheral B-cell populations in vivo. BLYS is made by immune-cells called monocytes and macrophages. When monocytes and macrophages are activated, BLYS is released and binds to a receptor found only on B cells.

**配方:** Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

**浓度:** 1 mg/mL.

**保存:** +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

**应用:** ELISA, WB, ICC/IF