



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

BMP7 Anti-Human Antibody (Clone 4E7)

单克隆抗体

产品货号: MA01119

第三版

**免疫原:**Recombinant human BMP7 (293-431aa) purified from E. coli

**亚型:**Mouse IgG1 heavy chain and kappa light chain

**克隆:**Anti-human BMP7 mAb, clone 4E7, is derived from hybridization of mouse SP2/O myeloma cells with spleen cells from BALB/c mice

**描述:**BMP-7 (bone morphogenetic protein-7), also known as OP-1 (osteogenic protein-1), is a member of the transforming growth factor (TGF ) superfamily. BMP-7 plays an essential role during ontogeny of the mammalian eye and kidney, also, is expressed in heart, lung, small intestine, limb bud and teeth.

**配方:**Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

**浓度:**1 mg/ml

**保存:** +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

**应用:**ELISA, WB, ICC/IF