



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

c-MPL/TPOR (AMM2) Anti-Mouse Rat IgG MoAb Biotin

单克隆抗体

产品货号: MA01130

第三版

免疫原: Recombinant protein of mouse c-MPL extracellular domain.

亚型: IgG1

克隆: AMM2

描述: This receptor is activated when thrombopoietin protein binds to it, and the activated receptor stimulates a signaling pathway called the JAK/STAT pathway, which transmits signals from outside the cell to the cell's nucleus and is important for controlling the production of blood cells. This antibody specifically reacts with mouse c-MPL/TPOR protein.

配方: Lyophilized product in PBS containing 1 % BSA and 0.05 % NaN₃

浓度: Please see datasheet.

保存: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

应用: ELISA, WB, ICC/IF