

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

C/EBP- Anti-Human Antibody (Clone 47A1)

单克隆抗体

产品货号: MA01132

第三版

免疫原: Recombinant human C/EBP- (1-271 aa) purified from E. coli

亚型: Mouse IgG1 heavy chain and kappa light chain

克隆: Anti-human C/EBP beta mAb, clone 47A1, is derived from hybridization of mouse SP2/O myeloma cells with spleen cells from BALB/c mice

描述: CCAAT/enhancer-binding proteins (C/EBPs) are a family of basic leucine zipper (bZIP) DNA binding proteins that play important roles in normal cellular differentiation and function in a variety of tissues. The C/EBP family is critical for innate immunity and inflammatory response, and consists of five core members: C/EBP β , C/EBP δ , C/EBP α , C/EBP γ , and C/EBP ϵ .

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1 mg/ml

保存: +4 °C 保存 (1-2 周). 长期保存在 -20 °C 或者 -70 °C. 避免反复冻融.

应用: ELISA, WB, ICC/IF