

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Casein Kinase 1 Anti-Human Antibody (Clone AT2E2)

单克隆抗体

产品货号: MA01139

第三版

免疫原:Recombinant human Casein Kinase 1 alpha (1-337aa) purified from E. coli

亚型:Mouse IgG2b heavy chain and kappa light chain

克隆:Anti-human Casein Kinase 1 alpha mAb, clone AT2E2, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice

描述:Casein Kinase 1 alpha belongs to the protein kinase superfamily, CK1 Ser/Thr protein kinase family, Casein kinase I subfamily. CK1 isoforms (alpha, beta, gamma, delta, epsilon) and their splice variants are involved in diverse cellular processes including membrane trafficking, circadian rhythm, cell cycle progression, chromosome segregation, apoptosis and cellular differentiation

配方:Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度:1 mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF