



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

CASQ2 Anti-Human Antibody (Clone AT4E10)

单克隆抗体

产品货号: MA01140

第三版

免疫原: Recombinant human CASQ2 (20-399aa) purified from E. coli

亚型: Mouse IgG1 heavy chain and kappa light chain

克隆: AT4E10

描述: CASQ2 (Calsequestrin 2) protein responsible for holding Calcium together with other two proteins junction and triadin at the sarcoplasmic reticulum renders the heart vulnerable to severe cardiac arrhythmias, known as CPVT (catecholaminergic polymorphic ventricular tachycardia), during physical exercise.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1 mg/ml

保存: +4 °C 保存 (1-2 周). 长期保存在 -20 °C 或者 -70 °C. 避免反复冻融.

应用: ELISA, WB, ICC/IF