
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

CAV1 Anti-Human Antibody (Clone AT4C1)

单克隆抗体

产品货号: MA01143

第三版

免疫原: Recombinant human CAV1 (1-104aa) purified from E. coli

亚型: Mouse IgG2b heavy chain and kappa light chain

克隆: AT4C1

描述: The identification of caveolin-1 (CAV1) as the main structural component of caveolae, together with the finding that CAV1 might serve as a molecular organizer for membrane multiprotein complexes involved in cellular traffic, endo- and transcytosis, cell adhesion and signal transduction prompted a new impulse in the research on these intracellular organelles.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1 mg/ml

保存: +4 °C 保存 (1-2 周). 长期保存在 -20 °C 或者 -70 °C. 避免反复冻融.

应用: ELISA, WB, ICC/IF