



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-CD134/OX40 (B7-B5) Mouse IgG MoAb

单克隆抗体

产品货号: MA01149

第三版

免疫原:Recombinant soluble OX40

亚型:IgG1

克隆:B7-B5

描述:CD134/OX40 is a receptor of OX40 ligand (OX40L) transactivated by Tax of HTLV-1. OX40 is expressed in activated CD4+ and CD8+ T cells. OX40 on the activated T cells have a role in proliferation of T cells and production of cytokines by binding with OX40L expressed on antigen presenting cells like dendritic cells or activated B cells. In addition, it also be involved in activated T cell invasion to inflamed area through the binding with OX40L.

配方:Lyophilized product from PBS containing 1 % BSA and 0.05 % NaN<sub>3</sub>

浓度:Please see datasheet.

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF