



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-CD25 (H-48) Mouse IgG MoAb

单克隆抗体

产品货号: MA01153

第三版

免疫原:Purified human IL-2 receptor

亚型:IgG1

克隆:H-48

描述:CD25 is IL-2 receptor  $\alpha$  chain (IL-2R $\alpha$ ) which is a low-affinity receptor of IL-2. CD25 forms high-affinity IL-2 receptors associating with CD122 ( $\beta$  chain) and CD132 ( $\gamma$  c chain). CD25 is expressed in activated T cells and B cells. Increase of CD25 in serum is observed in individuals suffer adult T cell leukemia or malignant lymphoma.

配方:Lyophilized product from PBS containing 1 % BSA and 0.05 % NaN<sub>3</sub>

浓度:Please see datasheet.

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF