

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Chromogranin A Anti-Human Antibody (Clone 3A5)

单克隆抗体

产品货号: MA01175

第三版

免疫原: Recombinant human Chromogranin A (19-131aa) purified from E. coli

亚型: Mouse IgG1 heavy chain and kappa light chain

克隆: Anti-human Chromogranin A mAb, clone 3A5, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice

描述: Chromogranin A (CGA) is the major member of the granin family of acidic secretory glycoproteins that are expressed in all endocrine and neuroendocrine cells. Chromogranin A regulates catecholamine storage and release through intracellular (vesiculogenic) and extracellular (catecholamine release – inhibitory) mechanisms.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

应用: ELISA, WB, ICC/IF