

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Cib1 Anti-Human Antibody (Clone 1D1)

单克隆抗体

产品货号: MA01177

第三版

免疫原: Recombinant human CIB1 (1-191aa) purified from E. coli

亚型: Mouse IgG2b heavy chain and kappa light chain

克隆: Anti-human CIB1 mAb, clone 1D1, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice

描述: CIB1 (Calcium and integrin binding 1) is regulatory protein with 50% homology to calmodulin and calcineurin B, that encodes a member of the calcium-binding protein family. The specific function of this protein has not yet been determined; however this protein is known to interact with DNA-dependent protein kinase and may play a role in kinase-phosphatase regulation of DNA end joining.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C 或者-70 ° C. 避免反复冻融.

应用: ELISA, WB, ICC/IF