



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

CIS3/SOCS-3 (19A5) Anti-Human Mouse IgG MoAb

单克隆抗体

产品货号: MA01179

第三版

免疫原:Synthetic peptide for a C-terminal part of Human CIS3/SOCS-3 (C204)

亚型:IgG1

克隆:19A5

描述:The Janus family of protein tyrosine kinases (JAKs) and STAT transcription factors regulate cellular processes involved in cell growth, differentiation, and transformation through their association with cytokine receptors. The CIS family of proteins (also referred as the SOCS or SSI family) has been implicated in the regulation of signal transduction by a variety of cytokines. The cytokine-inducible SH2 protein-3 (CIS3/SOCS-3/SSI-3)

配方:Lyophilized product from 1% BSA in PBS containing 0.05% NaN₃

浓度:200 ug/mL

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF