



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

COX-2 (13H14) Anti-Human Mouse IgG MoAb

单克隆抗体

产品货号: MA01191

第三版

免疫原:Synthetic peptide for a part of Human COX-2

亚型:IgG1

克隆:13H14

描述:Cyclooxygenase (COX) is a membrane bound enzyme responsible for the oxidation of arachidonic acid to Prostaglandin G2 (PGG2) and the subsequent reduction of PGG2 to PHG2.

配方:Lyophilized product from 1% BSA in PBS containing 0.05%NaN3

浓度:100 ug/mL

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF