



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

CRMP-2 (C4G) Anti-Human Mouse IgG MoAb

单克隆抗体

产品货号: MA01199

第三版

免疫原:Recombinant Human CRMP-2

亚型:IgG1

克隆:C4G

描述:CRMP-2 (Collapsin response mediator protein-2) is highly expressed in the developing nervous system and identified as protein involved in Semaphorin 3A (Collapsin-1) intracellular signaling.

配方:Lyophilized product from 1% BSA in PBS containing 0.05% NaN₃

浓度:100 ug/mL

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF