

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

DACT3 Anti-Human Antibody (Clone 2A5)

单克隆抗体

产品货号: MA01217

第三版

免疫原: Recombinant human DACT3 (89-202aa) purified from E. coli

亚型: Mouse IgG1 heavy chain and kappa light chain

克隆: Anti-human DACT3 mAb, clone 2A5, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice

描述: DACT3, a member of the DACT (Dpr/Frodo) gene family, plays an important role in a negative regulator of Wnt/-catenin signaling. Unlike other Wnt signaling inhibitors that are silenced by DNA methylation, DACT3 repression is associated with bivalent histone modifications.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1 mg/mL

保存: +4 °C 保存 (1-2 周). 长期保存在 -20 °C 或者 -70 °C. 避免反复冻融.

应用: ELISA, WB, ICC/IF