

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

DYNLL Antibody (Clone AT13C9)

单克隆抗体

产品货号: MA01224

第三版

**免疫原:** Recombinant human DYNLL (1-89aa) purified from E. coli

**亚型:** Mouse IgG2b, k

**克隆:** Anti-human DYNLL mAb, clone AT13C9, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human DYNLL protein.

**描述:** The LC8 family members of dynein light chains (DYNLL1 and DYNLL2 in vertebrates) are highly conserved ubiquitous eukaryotic homodimer proteins that interact, besides dynein and myosin 5a motor proteins, with a large (and still incomplete) number of proteins involved in diverse biological functions.

**配方:** Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

**浓度:** 1 mg/ml

**保存:** +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

**应用:** ELISA, WB, ICC/IF