



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

/E2 Junction (1A1) Anti-Human Mouse IgG MoAb

单克隆抗体

产品货号: MA01225

第三版

免疫原:Synthetic peptide of part of Human Tau/E2 Junction (KSTPTAEAEAGIG)

亚型:IgG1

克隆:1A1

描述:Tau protein is a microtubule-associated protein expressed in neuroaxon and has molecular weight of 50~70 kDa. Some Alzheimer's patients have sites where tau protein is accumulated abnormally

配方:Lyophilized product from 1% BSA in PBS containing 0.05%NaN<sub>3</sub>

浓度:100 ug/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF