



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

EIF2S1 Anti-Human Antibody (Clone AT5E10)

单克隆抗体

产品货号: MA01232

第三版

免疫原:Recombinant human EIF2S1 (1-315aa) purified from E. coli

亚型:IgG2b heavy chain and kappa light chain.

克隆:5E10.

描述:Eukaryotic translation initiation factor 2, subunit 1 alpha (EIF2S1) is plays a central role in initiating translation, and provides for a rate-limiting step in protein synthesis. The phosphorylation of alphasubunit of EIF2S1 effectively prevents formation of the EIF2S1/GTP/Met-tRNAⁱ complex and inhibits global protein synthesis.

配方:Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度:1 mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF