



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Endothelin-B Receptor (8Z11) Anti-Human Mouse IgG MoAb

单克隆抗体

产品货号: MA01234

第三版

免疫原:ET-1bound recombinant human 6hN-ETBR produced in insect sf-9 cells

亚型:IgG2a

克隆:8Z11

描述:ET-A receptor is ET-1 specific and expressed in lung, aorta, atrium, liver, muscle, cerebral cortex, cerebellum, kidney, adrenal, duodenum, colon and placenta at the mRNA or cDNA level, while ET-B receptor lacks isopeptide specificity and is found in liver, lung, placenta, kidney, heart and brain.

配方:Lyophilized product from 1% BSA in PBS containing 0.05%NaN3

浓度:100 ug/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF