



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

EphA2 Anti-Human Antibody (Clone 3D7)

单克隆抗体

产品货号: MA01238

第三版

**免疫原:** Recombinant human EphA2 (559-976aa) purified from E. coli

**亚型:** Mouse IgG2b heavy chain and kappa light chain

**克隆:** Anti-human EphA2 mAb, clone 3D7, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice

**描述:** The Eph family of receptors is the largest known subfamily of receptor tyrosine kinases (RTK). Eph RTKs and their ligands, the ephrins, are frequently overexpressed in a variety of cancers and tumor cell lines.

**配方:** Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

**浓度:** 1 mg/mL

**保存:** +4 °C 保存 (1-2 周). 长期保存在 -20 °C 或者 -70 °C. 避免反复冻融.

**应用:** ELISA, WB, ICC/IF