

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

EPM2A Anti-Human Antibody (Clone k2A3)

单克隆抗体

产品货号: MA01241

第三版

**免疫原:** Recombinant human EPM2A (243-331aa) purified from E. coli

**亚型:** Mouse IgG1 heavy chain and kappa light chain

**克隆:** Anti-human EPM2A mAb, clone k2A3, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice

**描述:** Epilepsy, progressive myoclonus type 2A (EPM2A), also known as laforin, is a dual-specificity phosphatase that associates with polyribosomes. The protein may be involved in the control of glycogen metabolism, particularly in monitoring for and preventing the formation of poorly branched glycogen molecules.

**配方:** Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

**浓度:** 1 mg/mL

**保存:** +4 °C 保存 (1-2 周). 长期保存在 -20 °C 或者 -70 °C. 避免反复冻融.

**应用:** ELISA, WB, ICC/IF