



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

ERK (Clone 5) Anti-Human Mouse IgG MoAb

单克隆抗体

产品货号: MA01242

第三版

免疫原:Recombinant Human ERK1

亚型:IgG1,

克隆:Clone 5

描述:The cell signaling molecules in Mitogen-Activated Protein (MAP) kinase family are playing the central role in the cell response. Extracellular Signal-Regulated Kinase (ERK), c-jun amino-terminal kinase (JNK), and p38 are well known molecules in the family.

配方:Lyophilized product from 1% BSA in PBS containing 0.05%NaN3

浓度:100 ug/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF