

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Fabp4 Anti-Human Antibody (Clone 3F4)

单克隆抗体

产品货号: MA01245

第三版

免疫原:Recombinant human FABP4 (1-132aa) purified from E. coli

亚型:Mouse IgG1 heavy chain and kappa light chain

克隆:Anti-human FABP4 mAb, clone 3F4, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice

描述:FABP4 (Fatty acid binding protein 4) is expressed in adipocytes and macrophages, and integrates inflammatory and metabolic responses. FABP4 is thought to regulate fatty acid uptake, release, and storage in adipocytes and participates in systemic glucose homeostasis and in macrophage responses in atherosclerosis.

配方:Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度:1 mg/mL

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF