

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Fabp7 Anti-Human Antibody (Clone At1D1)

单克隆抗体

产品货号: MA01246

第三版

免疫原: Recombinant human FABP7 (1-132aa) purified from E. coli

亚型: Mouse IgG2b heavy chain and kappa light chain

克隆: Anti-human FABP7 mAb, clone AT1D1, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice

描述: FABP7 (Fatty acid binding protein 7) is also known as brain fatty acid-binding protein (B-FABP) and is a member of fatty acid-binding proteins (FABPs) which are a family of small, highly conserved, cytoplasmic proteins to bind long-chain fatty acids and other hydrophobic ligands.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1 mg/mL

保存: +4 °C 保存 (1-2 周). 长期保存在 -20 °C 或者 -70 °C. 避免反复冻融.

应用: ELISA, WB, ICC/IF