

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

FKBP6 Antibody (Clone AT9B7)

单克隆抗体

产品货号: MA01255

第三版

免疫原: Recombinant human FKBP6 (1-327aa) purified from E. coli

亚型: IgG1,k

克隆: Anti-human FKBP6 mAb, clone AT9B7, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human FKBP6 protein.

描述: FKBP6 (FK506 binding protein 6) is a 327 amino acid protein. FKBP6 is a member of the immunophilin protein family, which play a role in immunoregulation and basic cellular processes involving protein folding and trafficking. FKBP6 is essential for homologous chromosome pairing in meiosis during spermatogenesis. FKBP6 is ubiquitously expressed in all tissues, with higher expressed in testis, heart, skeletal muscle, liver, and kidney. Deletion of FKBP6 may contribute to hypercalcemia and growth delay in Williams-Beuren syndrome.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1 mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

应用: ELISA, WB, ICC/IF