



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

G-CSF (4-12-2) Anti-Human Mouse IgG MoAb

单克隆抗体

产品货号: MA01262

第三版

免疫原:Recombinant Human G-CSF

亚型:IgG1

克隆:Apr 12, 0002

描述:G-CSF (Granulocyte Colony Stimulating Factor) is cytokine which accelerates differentiation and multiplication of neutrophic precursor cell.

配方:Lyophilized product from 1% BSA in PBS containing 0.05% NaN3

浓度:500 ug/mL

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF