



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Galectin-3 (87B5) Anti-Human Mouse IgG MoAb

单克隆抗体

产品货号: MA01268

第三版

免疫原: Recombinant Galectin-3 protein

亚型: IgG2a

克隆: 87B5

描述: Galectin is widely distributed in nematodes, insects, and porifers, as well as vertebrates, and it has also been found to be present in true fungi.

配方: Lyophilized product from 1% BSA in PBS containing 0.05% NaN₃

浓度: 100 ug/mL

保存: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

应用: ELISA, WB, ICC/IF