

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

GRO/CINC-1 (6G2) (C) Anti-Rat Mouse IgM MoAb

单克隆抗体

产品货号: MA01293

第三版

免疫原:Synthetic peptide to C-terminal of Rat GRO/CINC-1

亚型:IgM

克隆:6G2

描述:Cytokine-induced neutrophil chemo attractant 1(CINC-1) was originally purified from media conditioned by IL-1 stimulated rat kidney epithelioid cells (NRK-52E).

配方:Lyophilized product from 1% BSA in PBS containing 0.05% NaN₃

浓度:500 ug/mL

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF