

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

GSTT2 Antibody (Clone AT10H4)

单克隆抗体

产品货号: MA01301

第三版

免疫原: Recombinant human GSTT2 (1-244aa) purified from E. coli

亚型: Mouse IgG2b, k

克隆: Anti-human GSTT2 mAb, clone AT10H4, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human GSTT2 protein.

描述: Glutathione S-transferase theta-2 (GSTT2) is an enzyme catalyze the conjugation of reduced glutathione to a variety of electrophilic and hydrophobic compounds. Human GSTs can be divided into five main classes: Alpha, Mu, Pi, Theta, and Zeta. The theta class members GSTT1 and GSTT2 share 55% amino acid sequence identity and both are thought to have an important role in human carcinogenesis.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1 mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

应用: ELISA, WB, ICC/IF