



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

STEAP1 Anti-Human Antibody (Clone J2D2)

单克隆抗体

产品货号: MA01313

第三版

免疫原: Recombinant human STEAP1 (1-70aa) purified from E. coli

亚型: Mouse IgG2b heavy chain and kappa light chain

克隆: J2D2

描述: Six-transmembrane epithelial antigen of the prostate (STEAP) is expressed predominantly in human prostate tissue and in other common malignancies including prostate, bladder, colon, and ovarian carcinomas, and in Ewing's sarcoma, suggesting that it could function as an almost universal tumor antigen.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1 mg/ml

保存: +4 °C 保存 (1-2 周). 长期保存在 -20 °C 或者 -70 °C. 避免反复冻融.

应用: ELISA, WB, ICC/IF