



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Streptavidin Antibody (Clone 1C2)

单克隆抗体

产品货号: MA01315

第三版

免疫原:Recombinant Streptavidin (25-183aa) purified from E. coli

亚型:Mouse IgG2b heavy chain and kappa light chain

克隆:1C2

描述:Streptavidin, a tetrameric protein secreted by *Streptomyces avidinii*, binds tightly to a small growth factor biotin. It finds wide use in molecular biology through its extraordinarily strong affinity for the vitamin biotin.

配方:Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度:1 mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF