



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Synaptobrevin 2 Antibody (Clone 3E5)

单克隆抗体

产品货号: MA01316

第三版

免疫原:Recombinant human Synaptobrevin 2 (1-89aa) purified from E. coli

亚型:Mouse IgG1 heavy chain and kappa light chain

克隆:300000

描述:Synaptobrevin 2(Vesicle-associated membrane, VAMP2), which is an 18 kDa integral membrane protein localized to the cytoplasmic surface of synaptic vesicle, consists of a proline-rich N-terminal region, a highly conserved hydrophilic domain, followed by a transmembrane anchor and a C-terminal tail.

配方:Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度:1 mg/ml

保存: +4 °C 保存 (1-2 周). 长期保存在-20 °C或者-70 °C. 避免反复冻融.

应用:ELISA, WB, ICC/IF