



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

TARDBP Anti-Human Antibody (Clone k1B8)

单克隆抗体

产品货号: MA01317

第三版

**免疫原:** Recombinant human TARDBP (1-260aa) purified from E. coli

**亚型:** Mouse IgG1 heavy chain and kappa light chain

**克隆:** k1B8

**描述:** TAR DNA binding protein (TARDBP), also known as TDP-43, is a cellular protein. TARDBP was identified as the disease accumulating protein in patients with frontotemporal lobar degeneration (FTLD) with ubiquitin inclusions (FTLD-U) and in amyotrophic lateral sclerosis (ALS). TARDBP functions as a DNA-binding protein and specifically binds to the TAR DNA sequence motifs of HIV.

**配方:** Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

**浓度:** 1mg/ml

**保存:** +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

**应用:** ELISA, WB, ICC/IF