
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Tau (Phosphorylated) (C5) Anti-Human Mouse IgG MoAb

单克隆抗体

产品货号: MA01318

第三版

免疫原:PHF (paired helical filament) obtained from AD brain.

亚型:IgG1

克隆:C5

描述:Tau protein is a microtubule-associated protein expressed in neuroaxon and has molecular weight of 50 - 70 kDa. Some Alzheimer's patients have sites where tau protein is accumulated abnormally. This antibody reacts with phosphorylated part in amino acid sites, 386 - 406 (TDHGAEIVYKSPVVSVDTSR) of Tau protein. It does not react with the most carboxyl-terminal peptide (396 - 441) of normal human Tau protein.

配方:Lyophilized product from 1 % BSA in PBS containing 0.05 % NaN₃

浓度:

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF