

---

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

---

TBCA Anti-Human Antibody (Clone AT1A5)

单克隆抗体

产品货号: MA01320

---

第三版

**免疫原:**Recombinant TBCA (1-108aa) purified from E. coli

**亚型:**Mouse IgG2b heavy chain and kappa light chain

**克隆:**AT1A5

**描述:**Tubulin folding cofactor A, also known as TBCA, is one of four proteins (cofactors A, D, E, and C) involved in the pathway leading to correctly folded beta-tubulin from folding intermediates. Cofactors A and D are believed to play a role in capturing and stabilizing beta-tubulin in a quasi-native confirmation. This

**配方:**Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

**浓度:**1mg/ml

**保存:** +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

**应用:**ELISA, WB, ICC/IF