

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

TBCEL Antibody (Clone AT1B10)

单克隆抗体

产品货号: MA01322

第三版

免疫原: Recombinant human TBCEL (1-424aa) purified from E. coli

亚型: Mouse IgG2b heavy chain and kappa light chain.

克隆: Anti-human TBCEL mAb, clone AT1B10, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human TBCEL protein.

描述: Tubulin-specific chaperone cofactor E-like protein, also known as TBCEL, is a factor that regulates the microtubule cytoskeleton in determining cell behavior.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1 mg/ml

保存: +4 °C 保存 (1-2 周). 长期保存在 -20 °C 或者 -70 °C. 避免反复冻融.

应用: ELISA, WB, ICC/IF