

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

TEK/Tie-2 (2A1) Anti-Human Mouse IgG MoAb

单克隆抗体

产品货号: MA01325

第三版

免疫原: Recombinant protein of extracellular domain of human TEK/Tie-2

亚型: IgG1

克隆: 2A1

描述: TEK, also known as Tie-2 (Tyrosine kinase with Ig-like loops and Epidermal growth factor homology domains-2) and CD202b, is a 145 kDa receptor tyrosine kinase expressed on endothelial cells and hematopoietic cells.

配方: Lyophilized product from 1 % BSA in PBS containing 0.05 % NaN₃

浓度: Please see datasheet.

保存: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

应用: ELISA, WB, ICC/IF