

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Tenascin-C (Domain B) (4C8MS) Anti-Human Mouse IgG MoAb

单克隆抗体

产品货号: MA01327

第三版

免疫原: Human Recombinant Tenascin-C

亚型: IgG1

克隆: 4C8MS

描述: Tenascin-C (TN-C) is a component of the extracellular matrix (ECM) that has been shown to be involved in tissue interactions during fetal development and oncogenesis. It is glycoprotein consisting of six disulphide monomer isoforms being generated by alternative splicing.

配方: Lyophilized product from 1% BSA in PBS containing 0.05%NaN3

浓度: Please see datasheet.

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用: ELISA, WB, ICC/IF