

---

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

---

Tenascin-C (EGF Like Domain) (4F10TT) Anti-Human Mouse IgG

单克隆抗体

产品货号: MA01329

---

第三版

**免疫原:**Purified human glioma tenascin-C

**亚型:**IgG1

**克隆:**4F10TT

**描述:**Tenascin-C (TN-C) is a component of the extracellular matrix (ECM) that has been shown to be involved in tissue interactions during fetal development and oncogenesis. It is glycoprotein consisting of six disulphide monomer isoforms being generated by alternative splicing.

**配方:**Lyophilized product from 1% BSA in PBS containing 0.05%NaN3

**浓度:**Please see datasheet.

**保存:** +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

**应用:**ELISA, WB, ICC/IF