

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Tgfb1 Anti-Human Antibody (Clone K1H12)

单克隆抗体

产品货号: MA01331

第三版

免疫原: Recombinant human TGFBI (502-636aa) purified from E. coli

亚型: Mouse IgG2a heavy chain and kappa light chain

克隆: k1H12

描述: TGFBI, also known as BIGH3 and ig-h3 , is an extracellular matrix protein induced by transforming growth factor (TGF)-beta 1. TGFBI protein is involved in cell growth, cell differentiation, wound healing and cell adhesion.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1 mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

应用: ELISA, WB, ICC/IF